






EASIFLEX™ HA LIQUID SODIUM HYALURONATE

ACTIVE CONSTITUENT: 200mg/30mL HYALURONIC ACID (MW 1-1.25 million daltons)

Helps support joint function and your horse's mobility

EASIFLEX™ HA contains high molecular weight sodium hyaluronate which has anti-inflammatory properties, helps to maintain joint fluid elasticity and cartilage health, and aids in the lubrication of joints.

-  **HIGH QUALITY**
Each 15mL dose contains 100mg of top quality, high molecular weight Hyaluronic Acid
-  **LIQUID FORMULATION**
EASIFLEX™ HA is a liquid formulation which mixes well with horse feed
-  **APPLE/STRAWBERRY FLAVOUR**
Contains an apple/strawberry flavour which horses love
-  **EASY TO USE**
Both the 1L and 5L packs come with own dispensing cap
-  **PUMP PACK**
The 5L pack comes with a convenient 15mL pump pack



	EASIFLEX™ HA	Synovate HA®
Active Constituent	200mg/30mL HYALURONIC ACID	200mg/30mL HYALURONIC ACID
Molecular Weight	1-1.25 million daltons	1-1.25 million daltons
Route of Administration	Add to feed and mix well	May be added to feed or administered orally
Dose Rate / 500kg Horse		
Loading Dose (7-10 Days)	30mL / day	30mL / day
Maintenance Dose	15mL / day	15mL / day
Pack sizes	1L, 5L	946mL, 3.785L
Doses Per Pack at Maintenance Dose Rate	1L: 66 doses 5L: 333 doses	946mL: 63 doses 3.785L: 252 doses

EASIFLEX™ is a trademark of Abbey Laboratories Pty Ltd

Synovate HA® is a registered trademark of Kentucky Equine Research (Australasia) Pty Ltd

EASIFLEX™ HA LIQUID SODIUM HYALURONATE

ACTIVE CONSTITUENT: 200mg/30mL HYALURONIC ACID (MW 1-1.25 million daltons)



EASIFLEX™ HA contains high molecular weight sodium hyaluronate which has anti-inflammatory properties, helps to maintain joint fluid elasticity and cartilage health, and aids in the lubrication of joints.

FOR USE IN EQUINES ONLY

DIRECTIONS FOR USE

Maintenance: Administer 15mL per 500kg horse once daily.

Loading or pre-competition (optional): 30mL per 500kg horse may be given once daily for 7-10 days.

Add EASI-FLEX™ HA to feed and mix well.

CAUTION

Safe use in pregnant animals intended for breeding has not been proven.

WARNINGS

For animal use only.

Keep out of the reach of children and other animals.

In case of accidental overdose contact a health professional immediately.

This product should not be given to animals intended for human consumption.

STORAGE

Store below 30 °C (room temperature) in a cool, dry place out of direct light.

"I have been feeding EASIFLEX™ HA to my Eventers for the past 5 years. My 4-Star horse, Superclass is a testament to that, running consistently up the grades for 5 years to 4-Star level barely missing an event."

EASIFLEX™ HA is palatable and super easy to feed, my horses love it! The 5L size is very economical and comes with a handy pump dispenser which is very convenient given I am treating multiple horses daily."

MERRAN WALLIS – 5*Eventer

Learn more about Abbey Animal Health:



Find a store near you:



EASIFLEX™ is a trademark of Abbey Laboratories Pty Ltd